



Qubey-T

Fully Automatic Twin Scope Endowasher System





Minimize your Hygiene & Safety concerns. Disinfect your scope in the Qubey-T Endowasher after every procedure to maximize safety of the Patient, Staff and yourself.





eakage Tes,





Features

- Manual & Fully Automatic Operation Modes.
- Compatible with all Endoscopy systems.
- Special Tub design for minimum Disinfectant. USB Port for Data Transfer.
- Visible Disinfectant Level Indicator.
- Leakage Testing Facility.

- Spill Resistant User-Friendly Touch Screen Interface.
- In-built Printer for Reports.
- 50 Cycle Storage Memory.
- Powerful Channel Drying Mechanism.

Specifications - Ottomed Fully Automatic Endowasher (Qubey-T)

Ambient temperature		10°c – 40°c
Relative Humidity		30% to 85%
Water Supply Requirement	Flow Rate	12L/min – 16L/min
	Pressure	0.1MPa – 0.4 MPa
Main Unit Specification		
No. of Accommodated Endoscopes		2 nos.
Disinfectant Tank Capacity		15L
Leakage Testing		Bubble detecting during immersion (Visual)
Operating Modes		
Auto Mode		
Manual Mode	Operations	Time Range
	No. of wash cycle	1 – 20 cycles
	Disinfectant time	1 – 99 min
	Dry time	1 – 10 min
Alarm		At Low level
		End of process
Display		Touch Screen Panel
Memory Storage		50 cycles
Dimension		600 x 600 x 600mm (L x W x H)
Weight		55 Kg Approx
Power Supply		
Input Voltage		220 V
Allowable Voltage Fluctuation		±10%
Input Current		1.5 A
Frequency		50Hz
Allowable Frequency Fluctuation		±3









Manufactured by:

MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com